

Opal[®]

INSTRUMENT GRADE - HIGH LEVEL DISINFECTANT

OPAL[®] is an Instrument Grade - High Level Disinfectant solution for the reprocessing of heat sensitive semi-critical medical devices for which sterilization is not suitable. OPAL[®] is intended for use in manual reprocessing systems.

Please read Product Insert, MSDS and Technical Bulletin prior to using OPAL[®] Solution for information including warnings, contraindications and precautions.

STORAGE CONDITIONS & EXPIRATION DATE:

1. Store OPAL[®] in its original sealed container. Store below 30°C out of direct sunlight.
2. Once opened, the unused portion of the solution may be stored in its original container for up to 75 days until used.
3. The expiration date of OPAL[®] may be found on the container. DO NOT use product from an unopened or opened container after the labelled expiration date.
4. The reuse period of OPAL[®] should never exceed 28 days.

FIRST AID INFORMATION:

1. Blood, other body fluids and lubricants must be thoroughly cleaned from the surfaces and lumens of semi-critical medical devices before reprocessing in the disinfectant.
2. Rinsing of all the detergent, dry the probe fully.
3. Taking the dried probe, place it in the OPAL[®] solution.
4. Contact time for OPAL[®] is 6 min.
5. After 6 minutes, remove probe from OPAL[®] and rinse with a large volume of sterile water. Repeat this process three times. Always use fresh sterile water for each rinse. Alternatively, run the probe under a continuous stream of water for a minimum of 5 minutes.
6. Dispose of rinsing water.

Additional information is listed in the Material Safety Data Sheet available at www.whiteley.com.au.