

INTROCAN SAFETY® 3

CLOSED IV SYSTEM



Introcan Safety® 3 Product Specifications

	Size	Needle Length Inch	Needle Length mm	Flow Rate ml/min	Max Flow Rate at 300psi, ml/min*	Product Code
	24G	3/4"	19	22	4.5ml/sec	4251127-03
	22G	1"	25	35	7ml/sec	4251128-03
	20G	1"	25	65	12.5ml/sec	4251129-03
	20G	1 1/4"	32	60	12.5ml/sec	4251130-03
	18G	1 1/4"	32	105	17ml/sec	4251131-03
	18G	1 3/4"	45	100	17ml/sec	4251132-03

* Viscosity of Contrast Media = 10.4mPa

