

# Mayo Trolley

**TASK**  
MEDICAL

The Task Medical Mayo Trolley is made from stainless steel construction. Some of the key components of the mayo trolley are its removable tray, adjustable height and lockable castors.

## Height Adjustable

- Brush finished 304 stainless steel
- The top tray is detachable
- Four castors with brake
- Size 60 X 40 X 750 - 1150mm
- Tray Size 400 X 600mm
- Carrying capacity 35kg
- Dia 75mm castors

*Shallow tray  
25mm*

***Hospital Grade  
Stainless Steel***



[www.taskmedical.com.au](http://www.taskmedical.com.au)



SKU  
TKMAYOTR