

# GAMMEX<sup>®</sup> Latex Sensitive

S U R G I C A L   G L O V E S

## Latex, powder-free surgical gloves designed for superior tactile sensitivity and dexterity

- **Heightened tactile sensitivity:** The soft formulation combined with a unique process has resulted in a powder-free latex glove that is 20% thinner than Ansell's standard latex range offering an ideal balance of tactile sensitivity and dexterity
- **Elevated comfort:** These sterile surgical gloves also offer exceptional wearer comfort and protection thanks to its soft yet strong formulation
- **Enhanced grip:** The external surface is smooth with a micro-textured finish for secure handling of surgical instruments even in wet conditions
- **Reduced eye fatigue:** These latex gloves are brown in color to help prevent eye fatigue from surgical glare
- **Tested for chemo handling:** GAMMEX<sup>®</sup> Latex Sensitive powder-free gloves are tested in accordance with ASTM D6978 and can be safely used with a range of chemotherapy drugs\*\*

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

\*\*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

### Recommended For

- All surgeries, especially delicate procedures in ophthalmology, micro, and cardiovascular
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Extra comfort during long, sustained procedures
- Prevention of eye fatigue

### FEATURED TECHNOLOGIES



### GAMMEX<sup>®</sup> Innovation that Protects

Driven by innovation, GAMMEX<sup>®</sup> surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Latex Sensitive

SURGICAL GLOVES

## PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Brown
Thickness	Micro
Cuff Style	Straight
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Chlorinated
Allergy Prevention	None
Grip Level	Moderate
Double Gloving Recommendation	Outer glove or underglove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Diethyldithiocarbamate (ZDEC)

## ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	330051055	330051060	330051065	330051070	330051075	330051080	330051085	330051090

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com