


diagmed



diagnostic & medical

# ***Your Radiology Guide***

**One source, endless support**

*All your medical essentials, delivered when you need it.*

 **1300 788 928**

 [info@diagmed.com.au](mailto:info@diagmed.com.au)  [www.diagmed.com.au](http://www.diagmed.com.au)







# About Us

An Australian-owned SME and a leading Australian medical distributor, Diagmed is committed to enhancing healthcare through a vast product range and unwavering dedication. Specialising in practice set-ups and cost-effective solutions without compromising quality, we prioritise customer value and satisfaction with over 30 years of dedicated service and efficient distribution across Australia.

## Vision

To lead the way in advancing healthcare excellence throughout Australia by offering innovative solutions, superior products, and dedicated support, helping foster a healthier and more sustainable future.



### Delivery Excellence

*Commitment to consistently delivering superior products and services with the highest standards of quality, innovation, and reliability to ensure optimal healthcare outcomes.*



### Continuous Improvement

*Embrace innovation, actively adapt to the evolving healthcare landscape, and continually seek opportunities for improvement to stay ahead and enhance efficiency in healthcare delivery.*

 **1300 788 928**

 [info@diagmed.com.au](mailto:info@diagmed.com.au)  [www.diagmed.com.au](http://www.diagmed.com.au)



## Mission

We are dedicated to improving healthcare outcomes by delivering superior products, outstanding service, and dedicated support to the diverse landscape of Australian healthcare providers. Specialising in setting up new practices, we continuously adapt to meet evolving industry needs through our wide product range and cost-effective solutions. Our mission is to empower our customers with efficiency, savings, and quality, ensuring their success and contributing to overall healthcare improvement.

## Values



### Customer-Centric Approach

*Prioritise customer satisfaction through personalised solutions, dedicated support, and user-friendly platforms to enhance their experience and foster long-term relationships.*



### Integrity in Action



*Conducting business with honesty, transparency, and ethical principles, building trust and accountability among stakeholders, and upholding the integrity of the healthcare sector.*



# Contents

<b>Infusion Therapy</b>	<b>5</b>
• Biopsy Needles	6
• Spinal Needles	9
<b>Contrast</b>	<b>10</b>
• Contrast Agents	11
• Contrast Warmers	12
<b>CT Grids</b>	<b>13</b>
• CT Grids	14
<b>MRI</b>	<b>15</b>
• Furniture	16
• Consumables	17
<b>Radiation Protection &amp; Accessories</b>	<b>18</b>
• Apron & Shields	19
• Protective Eyewear	20
<b>Ultrasound</b>	<b>21</b>
• Ultrasound Gels	22
• Consumables	23
<b>Medical Imaging Accessories</b>	<b>24</b>
• Positioning Pads	25
<b>Medical Furniture</b>	<b>27</b>
• Trolleys	28
• Beds & Stools	30
• Stools	31
<b>Miscellaneous</b>	<b>32</b>
• Consumables	33
• X-Ray Stencil Markers	34
• Checklist new practice form	35

 **1300 788 928**


 [info@diagmed.com.au](mailto:info@diagmed.com.au)  [www.diagmed.com.au](http://www.diagmed.com.au)





diagmed

diagnostic & medical





# *Infusion Therapy*

**1** | Biopsy Needles  
Page 6-8

**2** | Spinal Needles  
Page 9

 **1300 788 928**

 [info@diagmed.com.au](mailto:info@diagmed.com.au)  [www.diagmed.com.au](http://www.diagmed.com.au)



## Max-Core™ Disposable Core Biopsy Instrument

The Max-Core™ Disposable Core Biopsy Instrument offers the convenience of a disposable with the simplicity of one-handed cocking, providing a flowing, consistent experience. All needles are engineered for exceptional performance, featuring ultra-sharp tips and polished surfaces for smooth entry and exit. A broad selection of biopsy systems and needle sizes assures an ideal match with procedural needs and personal preferences.

- Ergonomic handle for easy handling and control, with one-handed cocking
- 22 mm penetration depth for precise performance
- Two firing buttons for user preference
- Colour-coded design for needle gauge identification
- Compatible with TruGuide™ Coaxial Biopsy Needle for versatility

Description	Order Code
14G X 10CM	MDBRMC1410
16G X 10CM	MDBRMC1610
16G X 16CM	MDBRMC1616
18G X 10CM	MDBRMC1810
18G X 16CM	MDBRMC1816
18G X 20CM	MDBRMC1820
20G X 10CM	MDBRMC2010
20G X 20CM	MDBRMC2020



## Marquee™ Disposable Core Biopsy Instrument

The Marquee™ Disposable Core Biopsy Instrument offers two firing modes and two penetration depth options, providing you with versatile solutions tailored to address various soft tissue biopsy challenges.

- Unparalleled procedural versatility from a comprehensive core needle
- Powerful automatic fire mode
- Dual mode disposable core biopsy device with adjustable penetration depth (18 mm and 25 mm)
- Semi-automatic mode for firing the stylet only provides visual confirmation of sample notch placement

### Marquee™ Disposable Core Biopsy Instrument Kit With Coaxial Needle

Description	Order Code
14G X 13CM	MDBRMQ1413-KIT
16G X 10CM	MDBRMQ1610-KIT
18G X 10CM	MDBRMQ1810-KIT

Description	Order Code
14G X 10CM	MDBRMW1410
16G X 10CM	MDBRMQ1610
18G X 10CM	MDBRMQ1810
18G X 16CM	MDBRMQ1816
14G X 10CM	MDBRMQ1410






### **Mission™ Disposable Core Biopsy Instrument**

The Mission™ Disposable Core Biopsy Instrument seamlessly integrates semi-automatic functionality with intuitive visual status indicators and enhanced ultrasound visibility, ensuring effortless usage and precise tissue sampling.

- Compact design fits CT gantry
- Ergonomic grip for versatility
- Adjustable throw (10 mm and 20 mm) with depth indicator
- Visual Fire Ready Indicator ensures sample notch advances

Description	Order Code
16 X 10CM	MDBRMS1610
18G X 10CM	MDBRMS1810
20G X 10CM	MDBRMS2010
18G X 10CM	MDBRMSK1810
18G X 16CM	MDBRMSK1816
18G X 20CM	MDBRMSK1820
20G X 10CM	MDBR2010M
20G X 20CM	MDBRMSK2020




### **UltraCor™ Twirl™ Breast Tissue Marker**

The UltraCor™ Twirl™ Breast Tissue Marker's distinctive twirled ring shape and profile aids in its visibility.

- Distinctive twirled ring shape for enhanced visibility
- Polymer and gel-free construction for reduced introduction of materials into the breast
- Self-incorporating design for accurate placement directly into tissue
- Rigid, easy-to-use applicator with a 17G needle and bevelled tip for ultrasound visualisation
- Deploys independently or through a coaxial needle system

Description	Order Code
17G	MDUCTW17



### **UltraClip™ Dual Trigger Breast Tissue Marker**

The UltraClip™ Dual Trigger Breast Tissue Marker offers two feminine marker shapes, accommodating your patient with a choice of a front trigger and an alternate rear trigger to accommodate user preference, all while providing an ergonomic design for a comfortable, non-slip grip.

- Five PVA-enhanced shapes for ultrasound enhancement
- Includes Heart and Venus shapes for marking versatility
- MRI-friendly with options from three metals
- Exclusive dual triggers for easy deployment
- Accurate placement with 17G rigid needle design

Description	Order Code
17G X 10CM	MDUCLIP10
17G	MDUCLIP17
17G	MDUCLIP17DU
Ribbon Shape	MDUCLIPUS
Coil Shape	MDUCOIL17
17G Dual Trigger	MDUCOIL17D
17G Dual Trigger	MDUCOIL17DL





## Evolution Biopsy Needle

Temno Evolution® biopsy devices maintain the key features of the original Temno™ while offering enhanced procedural flexibility and performance. They come with or without a coaxial introducer needle.



**TEMNO™**

Description	Order Code
14G X 11CM	MDTT1411
14G X 15CM	MDTT1415
14G X 6CM	MDTT1406
16G X 6CM	MDTT1606
16G X 11CM	MDTT1611
16G X 15CM	MDTT1615
16G X 20CM	MDTT1620
18G X 6CM	MDTT1806
18G X 11CM	MDTT1811
18G X 15CM	MDTT1815
18G X 20CM	MDTT1820
20G X 6CM	MDTT2006
20G X 11CM	MDTT2011

Description	Order Code
20G X 15CM	MDTT2015
20G X 20CM	MDTT2020
22G X 11CM	MDTT2211

## Original Temno™ Biopsy Device

The Original Temno™ biopsy device is easy to use, offering simplicity, reliability and precision you can count on. It features a 20 mm sample notch and a lightweight, compact design. Available with or without a coaxial introducer needle.



**TEMNO™**

Description	Order Code
14G X 6CM	MDT1406
14G X 11CM	MDT1411
16G X 6CM	MDT1606
16G X 9CM	MDT1609
16G X 11CM	MDT1611
18G X 6CM	MDT1806
18G X 9CM	MDT1809
18G X 11CM	MDT1811
18G X 15CM	MDT1815
20G X 9CM	MDT2009
20G X 11CM	MDT2011
20G X 15CM	MDT2015

## Temno Universal Coaxial Introducer Needle

Universal Coaxial Introducer Needles simplify multi-site biopsies by providing a clear path. Compatible with Achieve® or Temno core needles and fine needle aspiration (FNA) needles.



**TEMNO™**

Description	Order Code
14G X 11CM	MDCTT1411
16G X 11CM	MDCTT1611
16G X 15CM	MDCTT1615
16G X 20CM	MDCTT1620
18G X 11CM	MDCTT1811
18G X 15CM	MDCTT1815
18G X 20CM	MDCTT1820
20G X 11CM	MDCTT2011
20G X 15CM	MDCTT2015

Description	Order Code
20G X 20CM	MDCTT2020

## Automatic Biopsy Device

The Achieve® programmable automatic biopsy system is the one solution for your soft tissue biopsy needs. This lightweight, automated system offers spring-loaded action for fast, accurate penetration and gives you the precise control and quality sampling capability you need when working with calcified or fibrous lesions.



**ACHIEVE™**

Description	Order Code
14G X 11CM	MDACHIEVE1411
16G X 9CM	MDACHIEVE1609
18G X 9CM	MDACHIEVE1809
18G X 11CM	MDACHIEVE1811
18G X 15CM	MDACHIEVE1815




### Spinocan With Quincke Bevel Needle

Experience enhanced safety and efficiency in spinal procedures with our specialised needle, featuring a sharp bevel for minimal puncture force and a fitted stylet to reduce the risk of tissue coring and potential epidermoid tumour development. The included prism in the hub further facilitates easy identification of cerebrospinal fluid (CSF) flashback.

- Sharp bevel design ensures low puncture force, optimising patient comfort during insertion
- Fitted stylet reduces the risk of tissue coring and potential epidermoid tumour formation, prioritising patient safety during the procedures


Description	Order Code
22G X 3.5"	BB4507908
25G X 3.5"	BB4505905
25G X 5"	BB4505913



### Quincke Spinal Needle

They are available in different sizes with lengths ranging from 90 mm to 178 mm and gauges from 18 G to 27 G. The key/slot arrangement of stylet and cannula hubs on the BDTM Quincke spinal needles facilitate proper needle bevel orientation.

Description	Order Code
18G X 3.5"	BD405248
18G X 6"	BD408360
19G X 3.5"	BD405250
20G X 3.5"	NDLBD2003.5
20G X 6"	BD405211
22G X 2.5" - 10	BD405244
22G X 2.5" - 25	NDLBD2709
22G X 3"	BD405255
22G X 3.5"	BD405256
22G X 5"	NDLBD2205
22G X 7"	NDLBD2207
25G X 2"	BD405245
25G X 3"	BD405246



Description	Order Code
25G X 3.5" - 10	BD405138
25G X 3.5" - 25	BD405257
27G X 3.5"	NDLBD2709

### Spinal Needle Quincke

This is the standard spinal needle with a conventional Quincke point, featuring excellent cutting performance and easy manipulation, ensuring successful procedures. It is available in a wide variety of sizes at reasonable prices to support different purposes and procedures.

Description	Order Code
18G X 90MM	NDLUNI1890
19G X 90MM	NDLUNI1990
20G X 70MM	NDLUNI2070
20G X 90MM	NDLUNI2090
21G X 70MM	NDLUNI2170
21G X 90MM	NDLUNI2190
22G X 70MM	NDLUNI2270
22G X 90MM	NDLUNI2290
23G X 70MM	NDLUNI2370



Description	Order Code
23G X 90MM	NDLUNI2390
25G X 90MM	NDLUNI2590





diagmed



diagnostic & medical

# Contrast

**1** | Contrast Agents  
Page 11


**2** | Contrast Warmers  
Page 12

 **1300 788 928**

 [info@diagmed.com.au](mailto:info@diagmed.com.au)  [www.diagmed.com.au](http://www.diagmed.com.au)



**Biliscopin 100ML**



Biliscopin is a contrast agent utilised for X-ray examinations focusing on the liver, bile ducts, gallbladder, and areas associated with bile flow. Its active ingredient, iotroxate meglumine, facilitates enhanced visualisation of these anatomical structures during the imaging process.

Order Code
BILISCPIN

**Gastrografin Oral 100ML**



Gastrografin is a contrast medium containing iodine used for X-rays or CT scans of the gastrointestinal tract, administered orally or as an enema, it enhances visibility in the digestive system, aiding in the clear depiction of specific areas for diagnostic purposes.

Order Code
GASTROGRAFIN

**Ioscan**



10ML



100ML

IOSCAN is an oral solution tailored for CAT-SCAN examinations of the gastrointestinal tract. Before ingestion, it should be diluted with water. Comprising sodium amidotrizoate, sorbitol, and essential additives, it serves as a preferred alternative to barium sulphate. It facilitates precise and vivid imaging for diagnostic purposes.

Description	Order Code
Ioscan 10ML	IOSCAN10
Ioscan 100ML	IOSCAN100

**Ultravist**

Ultravist is an injectable contrast medium, also known as a dye, that contains iodine. It's utilised to enhance visibility in specific areas of the body during X-ray examinations. This enhancement aids in clearly visualising anatomical structures such as the kidneys, bladder, heart, blood vessels, or spinal cord, facilitating more accurate diagnosis and investigation by doctors.

Description	Order Code
Ultravist 300, 1000ML - 10	ULTRAV300/1100
Ultravist 300, 20ML - 10	ULTRAV300/20
Ultravist 300, 50ML - 10	ULTRAV300/50
Ultravist 300, 75ML - 10	ULTRAV300/75
Ultravist 370, 100ML - 10	ULTRAV370/100
Ultravist 370, 50ML - 10	ULTRAV370/50



## Nuline MLI Refrigerator Incubator

The Nuline Refrigerator Incubator, designed for laboratories, pharmacies, and hospitals, ensures precise temperature control for your inventory. Equipped with a microprocessor-controlled heater and a digitally managed refrigeration system, it shields contents from fluctuations. Features like fan-forced circulation guarantee uniform temperature distribution. Security is paramount with standard door locks.

- Fan-forced circulation
- Voltage-free connection
- Door lock for security
- Door ajar alarms
- Interior LED lighting
- Automatic condensate water disposal
- Digital display with built-in alarms
- Hi / Lo temperature readout

59L



29L



29L



59L



- Gross capacity: 59L (also available in 29L)
- Temperature range: 2°C to 40°C
- Voltage: 220V

Description	Order Code
29L	NLMLI29
59L	NLMLI59



Ensure your inventory's safety and stability with the Nuline Refrigerator Incubator.



diagmed

diagnostic & medical





# CT Grids

1

CT Grids  
Page 14


 **1300 788 928**

 [info@diagmed.com.au](mailto:info@diagmed.com.au)  [www.diagmed.com.au](http://www.diagmed.com.au)

## Disposable CT Biopsy Grid


The Task Medical Disposable CT Biopsy Grid is a simple, disposable tool that ensures accuracy for CT biopsy and drainage procedures, saving both time and money. It guarantees precision regardless of pathology type or body part.

- 1cm markings for precise needle placement
- Flexible material conforms closely to patient's anatomy
- Disposable to prevent cross-contamination
- Tape strips for non-slip positioning
- Allows accurate marking with a felt-tip or marking pen
- Radiopaque grid with an identifiable dot-pattern on each CT image



### Order Code

TKGRID




## GuideLines CT Biopsy Grid 117 & 217

GuideLines® CT Biopsy Grids enhance the precision of needle entry during CT-guided biopsies, minimising patient discomfort and the need for repeat adjustments. Featuring self-adhesive properties, these grids ensure accurate localisation of needle entry sites, eliminating misregistration of skin marks and the need for painful re-sticks and rescans.

- Self-adhesive design optimises needle entry site localization
- Prevents misregistration of skin marks, reducing the need for repeat adjustments
- All-over adhesive backing conforms to patient anatomy, eliminating buckling
- GuideLines® and Mark-Through GuideLines® ensure first-stick accuracy
- Exclusive low-density non-metallic lines for clear imaging without obscuring artefacts
- Uniform spacing of grid lines ensures precision and minimises radiation exposure
- Extra spacing between specific grid lines aids orientation during procedures



217



117


Description	Order Code
CT Biopsy Grid 117	BEEK117
CT Biopsy Grid 217	BEEK217

Beekley Medical®  
when your diagnosis must be right™



diagmed

diagnostic & medical



# MRI

1

Furniture  
Page 16

2

Consumables  
Page 17



**1300 788 928**



[info@diagmed.com.au](mailto:info@diagmed.com.au)



[www.diagmed.com.au](http://www.diagmed.com.au)



## Easy Roller Wheelchair

EasyRoller is a revolutionary nonmetal wheelchair crafted entirely from environmentally friendly thermoplastics. Specifically designed as an MR SAFE patient transport device, it streamlines hospital and MR laboratory operations. With its unique construction, EasyRoller ensures patient comfort while facilitating the work of MRI technicians and hospital porters, allowing seamless movement through hospitals and MRI rooms without disrupting equipment.

- No metal parts
- No loose parts
- Recyclable parts
- Easy manoeuvring
- Wide and comfortable seat area
- High comfort and stability
- Approved for users of up to 150kg
- Aesthetic design
- Ergonomic form developed in a therapeutic environment
- Maintenance free



Made from  
recyclable  
PE and PA  
thermoplastics

100%  
non-metallic -  
non-magnetic



Order Code

WHEEASYRL



reddot

design award



Transport patients directly  
to the MRI Machine!



MADE IN  
GERMANY



EasyRoller



## MRI Conditional Fixed Height Stretcher

The MRI Safe Conditional Fixed Height Stretcher is expertly crafted from pure stainless steel, ensuring non-magnetic properties and thorough compatibility with SIEMENS 3.0 Tesla systems. Ideal for NMR and CT rooms, as well as veterinary hospitals, it is compatible with MRI systems from leading brands including SIEMENS, GE, PHILIPS, TOSHIBA, CANON, and HITACHI. The stretcher comes with essential certificates, including an MR compatibility report, CE (MDR 2017/745), EU free sales certificate, and Italy SFDA registration certificate.

- Stainless steel construction
- Drop down side rails
- 5" Lock Casters, 2 with brake, 2 without brake
- Adjustable headrest
- 2" pad included

**Order Code**

TKMRIST



## MRI Headset Covers

These disposable headphone covers offer an effective solution for scenarios requiring shared headphone usage. Crafted from white non-woven spun polypropylene, they provide a hygienic barrier against infectious bacteria and germs. With a latex-free elastic rim, they ensure compatibility with a wide range of large-size headphones.

- Latex-free elastic rim for compatibility with various large-size headphones
- Prevents the transfer of infectious bacteria and germs for hygiene maintenance

**Order Code**

HEA




## TASKSAFE Earplugs

TaskSafe Earplugs are specifically designed to offer superior hearing protection in environments with loud noises, supporting the auditory well-being of medical professionals. This makes them ideal for use in MRI suites, where powerful machines can generate noise levels exceeding 105dB (A).

- Individually Polybagged in Pairs
- Hearing Protector Class 5 Tested to AS/NZS 1270

**Order Code**

TKEPU



diagmed

diagnostic & medical





# Radiation Protection & Accessories

**1** | Aprons & Shields  
Page 19

**2** | Wall Racks  
Page 20

**3** | Protective Eyewear  
Page 20

 **1300 788 928**

 [info@diagmed.com.au](mailto:info@diagmed.com.au)  [www.diagmed.com.au](http://www.diagmed.com.au)



### Coat Apron Maxlite With Back Protection

Drawing on three decades of experience, UniRay creates radiation protection apparel through collaboration with top hospitals worldwide. Designed in Germany and crafted in India, our products feature cutting-edge German technology for superior quality. MaxLite™ redefines protective sheeting, offering a 20-30% weight reduction compared to traditional options and ensuring comfort and safety during x-ray surgeries and in Cath Labs.



### Patient Apron Litelead

Our patented LITELEAD™ material sets a new standard as the lightest fully leaded material available. Formulated with fine lead particles and an advanced binding system, LITELEAD™ ensures superior protection without the added weight associated with conventional lead oxide-based materials.

- Designed for lower body protection
- Meets IEC 61331-3:2014 & IEC 61331-1:2014 standards

Description	Order Code
30 X 30	UF1P3311NB
40 X 40	UF1P3322NB
50 X 50	UF1P3333NB



### Ovarian Shield Lite



The Uniray Ovarian Shield offers crucial protection for the female reproductive system during X-ray procedures. These shields are available in two sizes: child and adult, catering to different patient needs.

Description	Order Code
Child	UF733402BAM
Adult	UF733401PR

### Gonad Shield Lite



The Uniray Gonad Shield offers crucial protection for the male reproductive system during X-ray procedures. These shields are available in two sizes: child and adult, catering to different patient needs.

Description	Order Code
Child	UF723403BAM
Adult	UF723401NB

### Thyroid Shield



The Uniray Thyroid Shield is crafted to provide complete protection with a snug fit. Each shield can be easily attached to the apron using a convenient push-button mechanism.

Order Code
UF523300NB



## Thyroid CT Shield

This thyroid shield offers comprehensive protection for the neck and sternum against ionising radiation, ensuring safety during medical procedures.

- Adjustable fastener for customised fit
- Velcro closure for easy and secure fastening



### Order Code

KX1D3T100

## Wall Racks

Uniray offers two sturdy and durable wall racks for lead aprons, providing convenient storage solutions for medical facilities.

- Constructed from heavy-grade stainless steel for enhanced strength and durability
- Equipped with three or five swivel hanger arms for secure storage of lead aprons
- Swivel hanger arms offer full side-to-side movement, maximizing space-saving storage
- Stainless steel construction ensures longevity and ease of use
- Ideal for medical facilities seeking efficient and organised lead apron storage solutions

### 3 Arms

### Order Code

UFA08001



### 5 Arms

### Order Code

UFA09001




## Nike Rabid - 0.75MM

Stylish, lightweight frames designed to provide clarity and maximum eye coverage.

### Order Code

KX1D3T100



## Nike Brazen - 0.75MM

Stylish, universal fitting frames with an additional nose bridge for comfort.

### Order Code

UF830001P



## Front & Side - 0.75MM - 0.5MM

Protects from frontal and lateral radiation.

### Order Code

UF820001



## Max 30 - 0.75MM

Sporty look-comfortable padding on the temples and bridge of the nose.

### Order Code

UF840001



## Front - 0.75MM

Full frontal radiation protection zones and tight seals against the cheekbones and temples.

### Order Code

UF10001



## Fitover - 0.75MM

A universal solution for complete clarity and protection-fits over regular eyewear.

### Order Code


UF860001





diagmed

diagnostic & medical





# Ultrasound

**1** | Ultrasound Gels  
Page 22

**2** | Consumables  
Page 23

 **1300 788 928**

 [info@diagmed.com.au](mailto:info@diagmed.com.au)  [www.diagmed.com.au](http://www.diagmed.com.au)

Ultrasound  
**gel**  
water soluble  
hypoallergenic



## Clarity Ultrasound Transmission Gel

Clarity Ultrasound Gel is a top choice for medical ultrasound procedures, endorsed globally for its exceptional performance. It boasts an acoustically correct formula that ensures clear imaging while being gentle on sensitive skin. With high viscosity and vacuum technology, it protects probes and enhances comfort during procedures.

- Hypo-allergenic and non-irritating
- Water-soluble and easy to clean
- Free from formaldehyde, salt, and toxic substances
- Made from safe materials
- Odour-free and pleasant to use
- Does not damage probes, ensuring reliability
- Produced with vacuum treatment for quality assurance

Description	Order Code
Gel 250ML	CLUSG250
Gel 5L Clear	CLUSG5CL
Gel 5L Blue	CLUSG5BL

Now comes with a pump!



## Aquasonic Ultrasound Transmission Gel



Aquasonic 100 sets the global standard for medical ultrasound transmission gel, trusted universally for diagnostic and therapeutic procedures requiring viscosity.

Description	Order Code
Gel 250mL	AQUA250
Gel 5 L Clear	AQUACL
Gel 5 L Blue	AQUABL



## Scan Ultrasound Gel



This aqueous ultrasound gel is acoustically superior, ensuring optimal transmission of ultrasound waves across a broad frequency range. It's recommended by top ultrasound equipment manufacturers globally, and its bacteriostatic properties ensure hygiene without formaldehyde.

Order Code
SCANP



## Ultrasound Gel Warmers



Clarity Gel Warmer swiftly heats ultrasound gel with precise temperature control. It ensures safety with a thermal cut-off feature, offers durability for wall mounting or countertop use, and includes a two-year warranty.

Description	Order Code
Single	TKGELWAR
Double	TKGELWAR2



## Latex-Free Probe Covers

Clarity Latex-Free Probe Covers are a revolutionary solution tailored for endocavity imaging on ultrasound systems. Crafted from latex-free synthetic polyisoprene, these probe covers provide a seamless fit across diverse ultrasound transducers, ensuring utmost flexibility for healthcare professionals during examinations.

- Ultra smooth surface for enhanced manoeuvrability
- Single-walled construction for optimal sensitivity
- Translucent natural colour for clear visualisation
- Latex-free composition eliminates allergic reactions
- Non-lubricated design for versatile application

**Order Code**

CLPCVLF



## Latex Probe Covers

Clarity Lubricated Probe Covers are specifically tailored for gynaecological ultrasonography, particularly for transvaginal ultrasound procedures.

- Rolled probe covers fit various transvaginal transducers easily
- Each box has 100 foil-sealed covers for optimal hygiene
- Made to transducer cover specifications for a perfect fit
- Crafted from top-grade latex for reliability
- Meets ISO 4070 for quality assurance
- Single-use, non-sterile for easy disposal

**Description**

**Order Code**

Lubricated

TKPCL

Non-Lubricated

TKPC



## Sterile Freezer Bags

The Clarity Sterile Freezer Bags are designed for single-use convenience, providing a reliable solution for storing items in a sterile environment. Each bag measures 30 x 45cm, offering ample space for various storage needs. Additionally, they come individually packed for ease of use and hygiene maintenance.

- Promotes hygiene and minimises contamination with easy disposal after each use
- Each bag (30 x 45cm) offers ample storage space for various items
- Two latex-free bands ensure a snug closure, fitting items of different shapes and sizes

**Order Code**

CLFBW





diagmed

diagnostic & medical




# Medical Imaging Accessories

**1** | Positioning Pads & Sandbags  
Page 25

**2** | Positioning Pads  
Page 26

 **1300 788 928**


 [info@diagmed.com.au](mailto:info@diagmed.com.au)  [www.diagmed.com.au](http://www.diagmed.com.au)



### Task Medical Positioning Pad

The Task Ultrasound Trans Vaginal Wedge Patient Positioning Pad enhances female pelvic ultrasound procedures. Made from high-density foam with a hospital-grade vinyl cover, it ensures optimal patient positioning on non-tilting couches. Easy to clean with a zippered cover for convenience.

- Antibacterial
- Waterproof
- Extreme durability
- Clean with standard clinical wipes



### Body Block Purple

The Body Block Purple Patient Positioning Device offers versatile support for the head, arm, leg, and foot, and is specifically crafted for hygiene and ease of cleaning. It is an ideal fit for medical, beauty, and exercise settings.

- Body Block caters to medical needs from head to feet
- Polyurethane construction ensures easy cleaning and durability



### Sandbags

Made of durable vinyl with a herringbone handle for easy grip, these sandbags are filled with metal-free sand and stitched with non-conductive material for x-ray safety.

- These versatile sandbags are ideal for assisting with X-ray positioning while also serving as weight bags for traction




Description	Order Code
Small 1KG	SANS
Medium 3KG	SANM
Large 5KG	SANL


### Positioning Pad Set 7 Pieces

Dunlop's high-density radiolucent foam positioning pads are meticulously crafted in Australia to revolutionise medical imaging. Engineered to perfection, these pads are designed to ensure precise positioning, eliminating any chance of patient movement during diagnostic procedures. With Dunlop's expertise and commitment to excellence, you can trust us to elevate your diagnostic imaging experience to new heights.

- 2 x 15 Degrees
- 1 x Circle
- 1 x 45 Degrees
- 2 x Small 45 Degrees
- 1 x Step



**15 Degrees Wedge Positioning Pad**




**Order Code**

POS15



**45 Degrees Positioning Pad**




**Order Code**

POS45



**Small 45 Degrees Positioning Pad**




**Order Code**

POS45S



**Large 45 Degrees Positioning Pad**



**Order Code**

POS45L



**Circle Positioning Pad 20CM**




**Order Code**

POSCIRC



**Step Positioning Pad**



**Order Code**

POSSTEP





diagmed

diagnostic & medical



# Medical Furniture

1 | Trolleys  
Page 28-29

2 | Beds & Stools  
Page 30

3 | Stools  
Page 31

📞 1300 788 928

✉ [info@diagmed.com.au](mailto:info@diagmed.com.au) 🌐 [www.diagmed.com.au](http://www.diagmed.com.au)

### Emergency Trolley Premium

With a sturdy ABS table mounting and 304 stainless steel plate, it ensures enduring quality. The detachable drawers offer easy cleaning and replacement, while ample storage options provide organised access to supplies. Equipped with smooth, noise-free castors and brakes, it supports heavy-duty use with up to 50 KG drawer load and 200 KG overall load.

- Durable construction with ABS and stainless steel
- Detachable drawers for easy maintenance
- Ample storage for organised access



**Order Code**  
TKEMTAL5R



### Emergency Trolley

The Task Emergency Trolley has been designed for Australian healthcare facilities. Stronger construction and larger wheels make mobility highly responsive and effortless. With a range of specific features, this trolley offers everything you need in an emergency at your fingertips - fast!

- Central key locking
- Strong alloy frame
- Soft closing drawers
- External power board
- Washable ABS body



**Order Code**  
TKEMT5R



### Task Medical Stainless Steel Trolleys



Task's Medical instrument trolleys are crafted from Hospital Grade 304 Stainless Steel, showcasing a blend of durability and elegance. These trolleys are meticulously engineered to reduce the build-up of dirt and bacteria, facilitating easy cleaning and maintenance. With meticulous attention to detail, they are designed for maximum efficiency, ensuring a streamlined process in healthcare environments.

#### Rails, No Draws Single



Measurements	Order Code
50 X 50X 97CM	TKOD50N

#### Rails, No Draws Double



Measurements	Order Code
50 X 49X 80CM	TKOD80N

#### 1 Drawer Stainless Steel Trolley



Measurements	Order Code
49 X 49 X 90CM	TK1D50N



**2 Drawer Stainless Steel Trolley**



**TASK**  
MEDICAL

Measurements	Order Code
50 X 50 X 90CM	TK2D50N

**3 Drawer Stainless Steel Trolley**



**TASK**  
MEDICAL

Measurements	Order Code
50 X 50 X 90CM	TK3D50N

**4 Drawer Stainless Steel Trolley**



**TASK**  
MEDICAL

Measurements	Order Code
50 X 50 X 90CM	TK4D50N

**Side by Side 2 Drawer Stainless Steel Trolley**



Measurements	Order Code
80 X 50 X 90CM	TK2D80N

**TASK**  
MEDICAL

**Side by Side 4 Drawer Stainless Steel Trolley**



Measurements	Order Code
80 X 50 X 90CM	TK4D80N

**TASK**  
MEDICAL

### Task HI-LO Medical Examination Couch 3 Sections

This Examination Couch model is one of our most popular ones for Medical Practices. It is a no-compromise solution to work area accessibility, patient comfort, and operational flexibility and affords a full range of positioning options, from upright exam, to chair, to flat table.

- 3-Section 1- Motor Examination Couch
- Single electric actuator motor
- Height range - 480mm to 920mm
- Heavy-duty PVC vinyl that is fire, stain, mildew, oil, water and wear-resistant
- Retractable Castors and Levelling Feet
- 1980mm (Length) x 700mm (Width)
- Face Hole, Gas Strut Operated Sections
- Backrest tilts 90° to Flat (no down tilt)
- 250kg Safe Working Load
- Warranty: Frame - 5 Years, Electrical - 3 Years, Upholstery - 1 Year

Description	Order Code
Blue	TK0301NB
Black	TK0301BL
Grey	TK0301GY



### Winbex Stools

Winbex stools are designed for professionals in sedentary roles, prioritising support and freedom of movement. Dedicated to advancing ergonomics, particularly in the medical field, Winbex ensures optimal comfort and productivity for users.

#### Ascot Standard Saddle Seat



Designed and built for ergonomics and comfort, the backrest provides more support for longer working hours

**WINBEX**

Description	Order Code
Black	WX5105S02BL
Navy Blue	WX5105S02NB
With Back, Black	WX511114BL
With Back, Navy Blue	WX511114NB



#### Kimi Saddle Seat Stool



Helps to maintain correct posture, relieves stress from knees and hips, the backrest provides more support for longer working hours

**WINBEX**

Description	Order Code
Black	WX511201BL
Navy Blue	WX511201NB
With Back, Black	WX511202BL
With Back, Navy Blue	WX511202NB





### Doctor Stool T Seat

- Comfortable seat and backrest supports body and spine
- Adjustable levers for seat height and back tilt
- customisation
- Robust construction ensures durability for prolonged use



#### WINBEX

Description	Order Code
Black	WX510710BL

### Duma Doctor Stool

- Comfort seat and backrest supports body and spine
- Adjustable levers for seat height and back tilt customisation
- Robust construction ensures durability for prolonged use



#### WINBEX

Description	Order Code
Navy Blue	WXDUMANB

### Microscope Chair

- Separate adjustment for seat height, seat tilt and back angle
- Ergonomic seat design featuring an elegant curve for optimal comfort
- Equipped with smooth roller castors for effortless mobility



#### WINBEX

Description	Order Code
Black	WX5113301BL

### Mica Saddle Seat Stool

- Allows the stool to follow the user's body movement, enhancing comfort and reducing strain
- Ensures stability and durability while maintaining a sleek appearance



#### WINBEX

Description	Order Code
Black	WX5113201BL
Navy Blue	WX5113201NB

### Comfortline Surgeon Stool

- Heavy-duty caster: 65mm diameter
- Adjustable lift: 140mm / 200mm
- Hospital-grade vinyl padded seat for comfort



#### WINBEX

Description	Order Code
Black	WX002908BL
Navy Blue	WX002908NB

### Premium Surgeon Stool

- Facilitates easy adjustment and manoeuvring, enhancing user convenience and usability
- Single lever for height adjustment



#### WINBEX

Description	Order Code
Black	WX002327BL
Navy Blue	WX002327NB





# *Miscellaneous*

**1** | Consumables  
Page 33

 **1300 788 928**

 [info@diagmed.com.au](mailto:info@diagmed.com.au)  [www.diagmed.com.au](http://www.diagmed.com.au)



## Task Mixing Cannula

The Task Medical Mixing Cannula is a sterile, semi-rigid tube designed for precise liquid mixing in healthcare, laboratory, and research settings. With a tapered tip for superior manoeuvrability, it ensures contamination-free operations. Compliant with medical industry standards, its single-use design and EO sterilisation guarantee safety.

- Compatible with standard luer slip syringes, it simplifies mixing processes
- Each pack contains 100 cannulas, conveniently packaged in peel packs

### Order Code

TKMIX



## Bite Block Covers

Task Disposable OPG Bite Block Covers are designed to enhance hygiene and patient safety during dental imaging procedures. With 300 single-use covers per box, they offer a convenient and effective solution for maintaining cleanliness throughout each procedure.

- Designed for most OPG X-ray machines
- Reduces cross-contamination risk
- Universally sized for quick setup
- Made with high-quality materials for reliability
- Each box contains 300 single-use disposable covers

### Order Code

TKBBC



## Ear Rod Cone Covers

Our Ear Rod Cone Covers are essential accessories for cephalometric ear cones, designed to enhance hygiene and convenience during medical procedures. Crafted with quality materials, these disposable covers offer peace of mind for practitioners and patients alike.

- Easy-to-use covers for hassle-free application and disposal after each use
- Compatible with most cephalometric ear cones, ensuring versatility across medical settings

### Order Code


EARCC



## Chrome X-Ray Stencil

Our Chrome X-Ray Stencil is a versatile tool for radiographers, offering clear and precise markings for left (L) and right (R) orientations in medical imaging. Available in both 20mm and 25mm square sizes, these chromium-plated markers ensure visibility and durability, enhancing accuracy in X-ray interpretations.

- Chromium-plated construction ensures reliability and longevity for consistent performance



Description	Order Code
L Chrome X-Ray Stencil	ST20L
R Chrome X-Ray Stencil	ST20R


\*Please use the form on the next page to order Chrome X-Ray Stencils.


# X-Ray Stencil Markers


Diagnostic & Medical Pty Ltd - Fax: 1300 788 938  
Unit 4, 119 - 127 Wicks Road Macquarie Park NSW 2113

Phone: 1300 788 928 Email: info@diagmed.com.au Web: www.diagmed.com.au

<p>MARKL/RINIT20</p> 	Initials	Qty	Initials	Qty	Initials	Qty	Initials	Qty

<p>Plain "R" &amp; "L" Markers</p> 	Qty	Qty

<p>Mammography Stencil Markers</p> 	Initials	Qty


--

Without Initials Qty
----------------------

Numbers	Qty

A full range of custom made stencil markers & saddle markers available. Most stencils are made to order and can be made to your requirements. **Please contact us for more details on 1300 788 928.**

## Word Stencil Markers

Initials	Qty	Initials	Qty	Initials	Qty	Initials	Qty	*Any Stencil available, Please list below
1st Dose		Compression		Injection		Oblique		
After Food		Dorsal		Inspiration		PA		
Ant		Erect		Int Rotation		Vertical		
AP		Expiration		Lateral Left		Weight Bearing		
Before		Ext Rotation		LPO		Tomography		
Before Injection		Flexion		Mobile		Supine		
CC		Control		MLO		Semi Erect		

Name/Organisation:			
Address:			
		Post Code:	
Contact Name:		Phone:	
Order No:		Date:	



# Checklist For New Practices

<b>Contact Name:</b>		<b>Phone No:</b>	
<b>Address:</b>			
<b>Email:</b>			

<b>Products</b>	<b>Qty</b>	<b>Comments (i.e. size, colour, specifications etc.)</b>
Ultrasound Couch		
Examination Couch		
Gynaecology Couch		
Recovery Bay Stretcher		
MRI Safe Stretcher		
MRI Wheelchair		
Saddle Seats		
Stainless Steel Trolleys		
Emergency Trolley		
Treatment Trolley		
Linen Skips		
Drip Stand		
Step Stools		
Hand Wash Dispensers		
Paper Towel Dispensers		
Glove Dispensers		
Wheelchairs		
Viewing Boxes		
Curtains and Rails		
Examination Lights		
Contrast Warmers		
Lead Aprons & Accessories		
Defibrillator		
ECG Machines		
Gel Warmer		
Autoclave		
Vital Signs/Patient Monitor		
Vaccine Fridge		
Privacy Screens		
Doctor's Chairs		
Waiting Rooms Chairs		
Receptionist Chair		
Probe Disinfection System		
Soaking Cup Kit & Probe Rack		
Positioning Pads & TV Sponge		
Sand Bags		
Sharp Trolley		
Patslide		
Blood Pressure Monitor		
Pulse Oximeter		
Resus Kit		
Scales		
Consumables Products		
Blood Collection Chair		



**Cost-effectiveness**

**Customer**

**Commitment to Quality**

**Satisfaction**

**Efficiency**

**Teamwork**


**Reliability**

**Innovation**

**Inclusivity**

**Empathy**

**diagmed**



**Integrity**

**Honesty**

**Long-term**

**Partnerships**

**Trustworthiness**

**diagnostic & medical**

**Leadership**

**Excellence**

**Dedication**

**Trust**  
**Continuous Improvement**

**Transparency**

**User-friendly Platforms**

**Dedicated Support**

**Delivery Excellence**

