

May-June 2024

EOFY *Sales*

**TAX TIME
DEALS**

**MUST
END
JUNE
30**

Connect with Us

 **1300 788 928**

 info@diagmed.com.au  www.diagmed.com.au

Task Medical Examination Couches

The Task Medical Hi-Lo Examination Couch in Navy Blue is essential for healthcare professionals, offering 2-section and 3-section configurations for maximum efficiency and comfort. The 70cm wide surface and motorised height adjustment ensure seamless positioning. The 3-section model includes gas strut-operated sections and a face hole for diverse needs. Built with a robust frame, this durable couch is a long-term investment for any medical practice.

TASK
MEDICAL

2 Section Electric Hi-Lo Exam Couch



2 Sections | 1 Motor
Foot Controller

- SWL 250kg
- 1930 x 700mm
- Height adjustable 460 - 910mm



Blue
\$1170
★★★★★★★★
+ GST



Freight charges
may apply

3 Section Exam Couch



3 Sections | 1 Motor
Foot Controller

**Most Popular
for GP!**

- SWL 250kg
- 1980 x 700mm
- Height adjustable 480 - 920mm



Blue
\$1190
★★★★★★★★
+ GST

Black
\$1190
★★★★★★★★
+ GST

Grey
\$1190
★★★★★★★★
+ GST



Freight charges
may apply



Task Medical Examination Couches Black Edition

The Task Medical Hi-Lo Examination Couch in Black Edition is designed for efficiency and adaptability, available in 2-section and 3-section configurations. The 70cm wide surface and motorised height adjustment ensure optimal positioning. The 3-section model offers a detachable headrest for various procedures. With a robust frame and sleek black upholstery, this couch is a durable and stylish choice for any medical facility.

TASK
MEDICAL

2 Section Electric Hi-Lo Exam Couch



2 Sections | 1 Motor
Foot Controller

- Versatile
- Easy height adjustments
- 70cm surface
- Professional look



2 Section
\$1270
★★★★★★★★
+ GST



Freight charges
may apply

3 Section Exam Couch



3 Sections | 1 Motor
Foot Controller

- Seamless adjustments
- Long-lasting
- Easy to clean
- Lockable wheels



3 Section
\$1290
★★★★★★★★
+ GST



Freight charges
may apply

Enhance patient care with
our **Task Black Edition
Examination Couches**.
Designed for both caregivers
and patients, they offer
functionality with a sleek
appeal.

*All prices subject to 10% GST (unless marked as GST Free). Prices do not include freight or delivery insurance. Offer ends June 30th, 2024, or while stocks last.

No Drawer No Shelves

Single

49 x 49cm
x 90 cm



TK0D50N

\$360

★★★★★★★★
+ GST

Double

80 x 50cm
x 90 cm



TK0D80N

\$449

★★★★★★★★
+ GST



Freight charges
may apply

One Drawer



TK1D50N

\$409

★★★★★★★★
+ GST

Two Drawer



TK2D50N

\$509

★★★★★★★★
+ GST



Freight charges
may apply

Three Drawer



TK3D50N

\$790

★★★★★★★★
+ GST

Four Drawer



TK4D50N

\$890

★★★★★★★★
+ GST



Freight charges
may apply

Stainless Steel Medical Trolleys

Two Drawer Double



TK2D80N

\$690

★★★★★★★★
+ GST

Four Drawer Double



TK4D80N

\$990

★★★★★★★★
+ GST



Freight charges
may apply

Mayo Instrument Trolley



595 x 395 x 25mm tray

TKMAYOTR

\$290

★★★★★★★★
+ GST

TASK
MEDICAL

White Powder Coated Trolleys

One Drawer



TK1DPC

\$360

★★★★★★★★
+ GST

95mm front
castors

125mm
deep drawer

- One or Two drawers
- Bottom shelf
- Soft closing magnetic drawer
- Smooth running lockable castors
- Top and bottom side rails
- Outstanding finish
- Powder coated finish

*High grade trolley with
smooth castors and a
outstanding finish, it is
a great alternative to
Stainless Steel models*

Two Drawer



TK2DPC

\$400

★★★★★★★★
+ GST



Freight charges
may apply

*All prices subject to 10% GST (unless marked as GST Free). Prices do not include freight or delivery insurance. Offer ends June 30th, 2024, or while stocks last.

Doctor's Chairs

Ideal for Practitioner's spending many hours each day while seated in their consultation rooms.

Star
ECYS111H

\$240

★★★★★★★★★
+ GST

120kg
rating



Tristar
ECYS444

\$290

★★★★★★★★★
+ GST

120kg
rating



Miami II
ECYS113A

\$390

★★★★★★★★★
+ GST

150kg
rating



Hercules
ECYS50H

\$590

★★★★★★★★★
+ GST

175kg
rating



Freight charges
may apply

Five Drawer Emergency Trolley

- Central key locking
- Strong alloy frame
- Soft closing drawers
- External power board
- Washable ABS body

TKEMT5LBR

\$990

★★★★★★★★★
+ GST



Freight charges
may apply

**With a range of specific
features, this trolley offers
everything you need in
an emergency at your
fingertips - FAST!**

Privacy Screen

Task Medical mobile folding room divider offers quick deployment and mobility, serving as a versatile, hygienic solution for workspace optimization and patient privacy.

- Cost-saving layout rearrangement
- Quick deployment and mobility
- Temporary or permanent divider
- Compact, foldable modesty screen
- Smooth, easy-to-clean surface
- Removable wheels for cleaning
- Antibacterial with silver-ion technology
- Compact: 180cm x 185cm, 30cm panels

Blue

TKHSPSB

\$1299

★★★★★★★★★
+ GST

White

TKHSPSW

\$1299

★★★★★★★★★
+ GST



Freight charges
may apply

Blood Collection Chair

Task Medical Blood Collection Chair is extra wide, comfortable and supportive. The chair includes an armrest and a foldable footrest.

- Comfortable and supportive
- Integrated step stool
- Extra wide seat : 55cm
- Width including armrest: 91.5cm
- Arm rest width: 18.5cm
- Length including footrest: 84cm
- Seated capacity 200kg
- Dimensions: W 92.5cm D 85cm H 108cm

Grey

TKBLCHGY

\$549

★★★★★★★★★
+ GST

Black

TKBLCHBL

\$549

★★★★★★★★★
+ GST

Blue

TKBLCHNB

\$549

★★★★★★★★★
+ GST



Available on
Black, Grey
or Blue



Freight charges
may apply

*All prices subject to 10% GST (unless marked as GST Free). Prices do not include freight or delivery insurance. Offer ends June 30th, 2024, or while stocks last.

Winbex Doctor Stool T Seat

WX510710BL
\$479
★★★★★★★★★
+ GST

Freight charges may apply

Adjustable levers for seat height and back tilt

Black



Winbex Duma Doctor Stool

WXDUMANB
\$599
★★★★★★★★★
+ GST

Freight charges may apply

Robust construction ensures durability for prolonged use

Navy Blue



Winbex Comfortline Surgeon Stool

WX002908NB
\$259
★★★★★★★★★
+ GST

Freight charges may apply

Navy Blue

WX002908BL
\$259
★★★★★★★★★
+ GST

Black



Winbex Premium Surgeon Stool

WX002327BL
\$225
★★★★★★★★★
+ GST

Freight charges may apply

Black

WX002327NB
\$225
★★★★★★★★★
+ GST

Navy Blue



Stools

Winbex Ascot Standard Saddle Seat Stool

WX5105S02NB
\$469
★★★★★★★★★
+ GST

Black

WX511114BL
\$509
★★★★★★★★★
+ GST

Navy Blue

Freight charges may apply

Black With Backrest

Navy Blue With Backrest

WX5105S02BL
\$469
★★★★★★★★★
+ GST

WX511114NB
\$509
★★★★★★★★★
+ GST



Winbex Kimi Saddle Seat Stool

WX511201NB
\$469
★★★★★★★★★
+ GST

Black

WX511202BL
\$509
★★★★★★★★★
+ GST

Navy Blue

Freight charges may apply

Black With Backrest

Navy Blue With Backrest

WX511201BL
\$469
★★★★★★★★★
+ GST

WX511202NB
\$509
★★★★★★★★★
+ GST



Winbex Mica Saddle Seat Stool

Black
WX51132001BL
\$449
★★★★★★★★★
+ GST

Navy Blue
WX5113201NB
\$449
★★★★★★★★★
+ GST

Safe Working Load 140kgs

- Natural sitting position
- Reduces stress on knees and hips
- Maintains correct posture
- Convenient height and position adjustments
- Easy tilt adjustments for customizable comfort

WINBEX

Freight charges may apply



*All prices subject to 10% GST (unless marked as GST Free). Prices do not include freight or delivery insurance. Offer ends June 30th, 2024, or while stocks last.

Easy Roller

Easy Roller is a revolutionary non-metal wheelchair crafted entirely from environmentally friendly thermoplastics. Specifically designed as an MR SAFE patient transport device, it streamlines hospital and MR laboratory operations. With its unique construction, Easy Roller ensures patient comfort while facilitating the work of MRI technicians and hospital porters, allowing seamless movement through hospitals and MRI rooms without disrupting equipment.

- No metal parts
- No loose parts
- Recyclable parts
- Easy manoeuvring
- Wide and comfortable seat area
- High comfort and stability
- Approved for users of up to 150kg
- Aesthetic design
- Ergonomic form developed in a therapeutic environment
- Maintenance free



red dot
design award



WHEEASYRL
\$2795
★★★★★★★★★
+ GST

Made from
recyclable PE and PA
thermoplastics



Freight charges
may apply



Transport patients directly to the MRI Machine!


Task Patient Transporter Bed Manual Hi-LO

The Task Medical EzyM Multipurpose Mobility Bed is an advanced manual wind hi-lo patient bed designed to meet the needs of any medical practice. Built primarily from high-impact ABS polymer and a high-strength alloy frame, the EzyM offers a significant upgrade from traditional steel-framed gurneys. It functions as both a patient transport stretcher and a hi-lo bed, making it versatile for use in medical practices, aged care facilities, and hospitals.

Use as:

- Examination couch
- Procedure table
- Emergency Stretcher
- Patient transporter
- Recovery room bed

Includes new **extra thick**
mattress upgrade



TKEYZM
\$1499
★★★★★★★★★
+ GST



Freight charges
may apply



86cm
56cm



Multipurpose mobility bed ideal for any
medical practice

*All prices subject to 10% GST (unless marked as GST Free). Prices do not include freight or delivery insurance. Offer ends June 30th, 2024, or while stocks last.

Vital Signs Monitors

Micro Spirometer



Vitalograph®

The updated Micro is a hand-held spirometry solution that delivers versatile, full spirometry testing in a range of environments

VG100053

\$1590

★★★★★★★★
+ GST

Includes

- Vitalograph Device
- Studio software
- Protective carry pouch
- USB Cable

*Mouthpieces required sold separately



Welch Allyn Connex Spot Monitor 71WE-6



Includes 1 adult cuff

Baxter
Welch Allyn

- Pro6000 Thermoscan ear thermometry
- Nonin SpO2 technology
- SureBP records in 15 sec'
- BP Averaging for hypertension detection
- EMR Connectivity
- 7" Colour Touch-screen

WA71WE-6

\$2390

★★★★★★★★
+ GST

*Requires a stand / not included



Spot Vital Signs Monitor Model 4400

- SureBP records BP in 15 sec'
- BP Averaging for hypertension detection
- SureTemp Plus oral axillary temp'
- Nonin SpO2, EMR Connectivity
- 7" Colour Touch-screen

WA44WT-6

\$2090

★★★★★★★★
+ GST

Baxter
Welch Allyn

Includes 1 adult cuff



*Requires a stand / not included

VSM Stands

71WE-6 Desk Stand

WA7000-DST

\$325

★★★★★★★★
+ GST

71WE-6 or 4400 Mobile Stand & Basket

WA7000-MS3

\$590

★★★★★★★★
+ GST

4400 Desk Stand

WA4400-DST

\$325

★★★★★★★★
+ GST

Baxter
Welch Allyn



*All prices subject to 10% GST (unless marked as GST Free). Prices do not include freight or delivery insurance. Offer ends June 30th, 2024, or while stocks last.

OMRON BP Monitors

HEM 7144T1 OMRON

- Irregular heartbeat detection
- Hypertension indicator
- Wide range cuff 22-42cm for medium to large arms



OMRON7144

\$109

★★★★★★★★★

+ GST



HEM 736T1 OMRON

- Measures and detects High Blood Pressure and AFIB at the same time
- Multi cuff range from 22cm to 42cm



OMRON736T1

\$180

★★★★★★★★★

+ GST



E-1201 ECG Machine

- Simultaneous 12 leads ECG acquisition & display
- Simultaneous 3 leads ECG printing with 80mm thermal paper



Pacemaker detection & built-in rechargeable battery

ECGMAC-EM1201

\$2595

★★★★★★★★★

+ GST FREE



EOFY Bonus

Riester Cardiophon 2.0 Stethoscope
Complimentary! rrp \$290

HEM 907 Digital BP Monitor

- Simple push button operation
- Accurate Intellisense technology
- Automatic cuff inflation and deflation

MADE IN JAPAN



OMRON107

\$720

★★★★★★★★★

+ GST



Includes 3 CUFFS
Small
Medium
Large

ProBP 2000 Digital BP Device

- One-button touch triggers the device to start capturing on cuff inflation
- Includes one Adult #11 FlexiPort cuff
- Meets AAMI standards for blood pressure accuracy



WAPBP2000

\$229 ea

★★★★★★★★★

+ GST

Baxter
Welch Allyn

DS66 Trigger Aneroid Sphyg'

Baxter
Welch Allyn

Robust design with Gear-free DuraShock™ technology to remain in calibration longer

- Laser-engraved dial face for accuracy
- 15 Year calibration warranty
- Withstands a 30" drop

WA5098-30

\$249

★★★★★★★★★

+ GST

With multi cuff kit



CP150 12 Lead Interpretive ECG

- Colour 7" touch screen
- ECG preview on screen
- One touch operation
- Thermal printer

Baxter
Welch Allyn

WACP150A

\$3295

★★★★★★★★★

+ GST FREE



Base model
CP150A-1EN66

*All prices subject to 10% GST (unless marked as GST Free). Prices do not include freight or delivery insurance. Offer ends June 30th, 2024, or while stocks last.

24 HR Ambulatory BP Monitor

With IHB & AFib Screening!

- All-in-one monitoring solution: Ambulatory (ABPM), Home (HBPM), Office (OBP), Automatic Night Blood Pressure (ANBP), Automatic Self Blood Pressure (ASBP) measurement/monitoring
- Adult cuff and carry pouch Bluetooth, USB output to PC or printer
- "AND" analysis software

ANDTM2441

\$1490

★★★★★★★★★
+ GST



Front LCD &
Top OLED
Display!

The TM-2440 monitor is a compact and lightweight ABPM weighing just 122 grams with proven accuracy ± 3 mmHg

- IHB & AFib Screening
- Automatically adjusts proper Cuff Pressure
- Bluetooth (BLE0 or USB cable) data output to computer or printer

ANDTM2440

\$1290

★★★★★★★★★
+ GST

AND
A&D Medical

Includes Software,
Adult Cuff & Carry
Pouch



OLED Display



Cardioline TouchECG HD

The ECG system is designed for Windows touch screen, PC's & laptops

- A fully diagnostic 12 lead ECG acquisition unit
- Wireless Bluetooth connection to Windows device. A wired option is also available
- Glasgow University Interpretation as standard
- 3 Year Manufacturers Warranty

CARDIOLINE



Software
included

ECG103550

\$1990

★★★★★★★★★
+ GST

Wireless
12 Lead
PC ECG

Clements Ceevac Portable Suction Pump

A versatile and reliable suction source designed for high vacuum performance and minimal maintenance.

- Runs on 240V mains, internal battery, or 12VDC automotive lead
- Max vacuum of -80kPa (-600mmHg) and 21 litres per minute flow rate
- Lightweight, compact, with 60-minute battery life: perfect for emergencies, home and wards

CLEMENTS



Vibration-free
operation

SUC81030

\$1195

★★★★★★★★★
+ GST

Aeon A100D 2Mhz Foetal Doppler

A high performance Foetal Heart Beat Detector with color OLED displaying FHR value

DOPPAEON

\$150

★★★★★★★★★
+ GST

AEON
TECHNOLOGY



- Real Time FHR & (3mins) FHR curve modes
- FHR Range: 60-120 bpm
- HR Range for Manual Count mode: 30-250 bpm
- FHR Accuracy: ± 1 bpm

*All prices subject to 10% GST (unless marked as GST Free). Prices do not include freight or delivery insurance. Offer ends June 30th, 2024, or while stocks last.

Welch Allyn Diagnostic Tools

PanOptic **PLUS** Ophthalmoscope & MacroView **PLUS** Otoscope

With iExaminer switch from optical to digital

WA71-PM3LXE

\$1490

+ GST

iExaminer
Smartphone
Bracket



- NEW PanOptic Eye Cup
- Li-ion **Plus** USB C Power Handle & Cable

PanOptic **Basic** Ophthalmoscope & MacroView **Basic** Otoscope

WA71-PM2LXU

\$1190

+ GST

- Li-ion **Basic** Power Handle
- USB Charging Cup and USB Cable

Baxter
Welch Allyn

Coaxial Ophthalmoscope & MacroView Basic Otoscope



- Universal Desk Charger
- Two Li-ion **Basic** Power Handles
- 240V power adaptor/lead

WA71-SM2LDX-6



\$1290

+ GST

PanOptic Ophthalmoscope & MacroView Basic Otoscope



- Universal Desk Charger
- Two Li-ion **Basic** Power Handles
- 240V power adaptor/lead

EA71-PM2LDX-6



\$1490

+ GST

Pro6000 Thermoscan Ear Thermometer

Designed to take quick & accurate readings

WA06000-200

\$270

+ GST

Baxter
Welch Allyn

Includes 1 small cradle



Ear Thermometer

WELCARE
EAR THERMOMETER



- Fast & accurate reading in one second
- Ear & object temperature detection
- Probe cover detection & eject function



Normal Elevated High

WCWET100

\$49 ea

+ GST

Infrared Forehead Thermometer

- Fast reading & high accuracy
- No-Contact thermometer
- Memory recall and °F/°C switchable
- Large screen with fever alarm

THERMJT

\$39

+ GST

SEJOY



*All prices subject to 10% GST (unless marked as GST Free). Prices do not include freight or delivery insurance. Offer ends June 30th, 2024, or while stocks last.

Seca Scales & Measures

Seca Electronic Digital Flat Scales

Seca 803



150kg

SECA803
\$118
★★★★★★★★★
+ GST

Seca 813



200kg
Extra wide

SECA813
\$190
★★★★★★★★★
+ GST

- Slimline design
- 100g graduations
- Rubber coated platform

seca
Precision for health

Samaritan Semi Automatic Defibrillators

Semi-automatic shock delivery with easy to follow visual and verbal guides

HEARTSN350P

\$1690
★★★★★★★★★
+ GST FREE

HEARTSN500P

\$2190
★★★★★★★★★
+ GST FREE



HeartSine

- Metronome for correct speed CPR
- 8 year warranty
- Shock proof case dust cover

Walmay Liquid Nitrogen Cryospray

- 0.35 Litre flask capacity
- 4 cone shaped spray tips for accurate visual aiming
- 1 curved spray tip



CRYOSPRAY
\$690
★★★★★★★★★
+ GST FREE

Walmay

Seca 354 Baby Scale

Electronic baby scales with fine graduation, also usable as flat scales for children

- Capacity 20kg
- Graduation 5g < 10kg > 10g
- Breast milk in take function

SECA354
\$249
★★★★★★★★★
+ GST



ChoiceMMed Pulse Oximeter
Complimentary! rrp \$69



Seca 769 electronic Column Scales with BMI Function & Measuring Rod

BMI with the press of just one key!

- Non slip surface on platform
- Easy roll transport castors
- Intelligent functions
- 100g graduations / lbs / kgs

SECA769
\$890
★★★★★★★★★
+ GST



seca
Precision for health




Skin Cancer Diagnosis Illuco IDS 1100 Dermatoscope

ILIDS1100
\$1045
★★★★★★★★★
+ GST



- Switchable polarized & non-polarized illumination
- 10X magnification and 22mm field of view
- 3 levels of brightness control



Malignant Melanoma

Basal Cell Carcinoma BBC

Squamous Cell Carcinoma SCC

Actinic Keratosis

ILLUCO

EOY Bonus




rrp \$89
8mm Contact Plate
Complimentary!

Includes universal smartphone adaptor & leather pouch with belt clip

*All prices subject to 10% GST (unless marked as GST Free). Prices do not include freight or delivery insurance. Offer ends June 30th, 2024, or while stocks last.

STAY CONNECTED



Unit 4, 119 - 127 Wicks Road,
Macquarie Park, NSW, 2113

1300 788 928

www.diagmed.com.au

info@diagmed.com.au

in Diagnostic & Medical Pty Ltd

***For the latest updates and exclusive deals, follow us
on LinkedIn and visit diagmed.com.au today!***

Stay informed with the newest advancements and special offers.
Join our community and never miss an update!