

Blood Collection Chair



Task Medical Blood Collection Chair is extra wide, comfortable and supportive. The chair includes arm rests and a foldable foot rest.

Extra Wide

- Comfortable & supportive with a high density soft foam cushion
- Integrated step stool
- Extra wide seat: 55cm
- Width including arm rest: 91.5cm
- Arm rest width: 18.5cm
- Length including foot rest: 84cm



Seated Capacity 200kg



Foldable foot rest

Available in Black, Blue & Grey



www.taskmedical.com.au



SKU
TKBLDCHBL
TKBLDCHNB
TKBLDCHGY